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Produktbeschreibung

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Instructions for use

Atrauman dressings are single-use, sterile, non-medicated, impregnated dressings for atraumatic wound treatment.

Composition

Support fabric made of hydrophobic polyester tulle, impregnated with a non-medicated ointment based on triglycerides (neutral fats). The ointment contains (INCI-assigned names): Caprylic/Capric/Miristic/Stearic Triglycerides; Bis-Diglyceryl Polycaprylate-P2.

Properties and mode of action

The soft, thin support material made of polyester allows close contact with the whole surface of the wound which is necessary for the drainage of wound exudate. The water-repellent polyester

filaments and the smooth surface structure of the tulle counteract adhesion to the wound. Atrauman allows for atraumatic and painless dressing changes. Atrauman's non-adhering properties are further enhanced by the ointment mass with which the product is impregnated. It has no occlusive effects and avoids causing maceration of the underlying skin. Excessive amounts of exudate are drained away into an absorbent secondary dressing placed on top of the Atrauman dressing.

By keeping the edges of the wound soft and supple, the non-medicated, non-sensitizing, paraffin-free ointment base of Atrauman, allows undisturbed epithelialization. Because it leaves no residues on removal, it conditions the wound for possible surgical procedures and facilitates wound cleansing.

Intended purpose
Atrauman are single-use, sterile, non-medicated ointment dressings suitable for the treatment of superficial acute and chronic wounds of any type. As a non-adherent wound contact layer, it is particularly suitable for preventing the secondary dressing from sticking to the wound bed and for keeping wound edges and surrounding skin supple. Due to the neutral properties of the non-medicated ointment mass, Atrauman is especially useful in dermatology as well as for patients with sensitive skin and those sensitive to certain medications.

Atrauman is suitable for the treatment on the human skin by professional users.

Mode of application
Take the dressing out of the peel pack with both product covers in place and, if needed, cut it with sterile scissors to the size of the wound. After removal of one of the product covers, place this side of the dressing onto the wound and remove the second cover. Place a sterile, absorbent dressing pad over the Atrauman dressing to absorb exudate. Unless otherwise prescribed by the physician or another healthcare professional, apply a new Atrauman dressing at each dressing change.

Atrauman can be left on the wound for up to 7 days. The wear time of the dressing depends on the condition of the wound and the level of exudate.

Note: If there are therapeutic reasons for not changing the dressing, the ointment may be absorbed when the reutilization period prevents nursing it a leur intégrité et à leurs performances. Informations sur demande.

Contraindications

Do not use Atrauman on patients who may be allergic to any of its ingredients.

Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Special precautions
In the absence of available data supporting the use of this dressing on sensitive population groups such as infants, children, pregnant or nursing women, and in the absence of data to the contrary, on these population groups this dressing should be used with caution following a clinician's recommendation.

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Before opening the product packaging, inspect it carefully for signs of damage. Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or nonuniform, torn or incomplete seals.

Incident reporting
For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Date of revision of the text: 2023-03-18
AT – PAUL HARTMANN Pty. Ltd. - Macquarie Park, NSW 2113
GB – PAUL HARTMANN Ltd. - Heywood/Greater Manchester OL10 2TT
ZA – HARTMANN South Africa - Northridge, 2169

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