HARTMANN

Veroval[®]

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Shenzhen AOJ Medical Technology Co., Ltd.
Room 301&4-5F, Block A, Building A, Jingfa Intelligent
Manufacturing Park, Xiaweiyyuan,
Gushu Community, Xixiang Street,
Bao'an District, 518126 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA
Email: info@aojmedical.com
Website: <https://www.aojmedical.com>
Tel: 86-755-2778 6026

EU REP

Share Info GmbH
Am Schulzentrum 12, 41564 Kaarst, Germany
Telephone: 0049 0179 5666 508
Email: EU-Rep@share-info.com

PAUL HARTMANN AG
Paul-Hartmann-Straße 12
89522 HEIDENHEIM, GERMANY
www.hartmann.info

CE
0 1 2 3

Dear customer,

We are pleased that you have chosen to purchase a pulse oximeter from HARTMANN. Veroval® Pulse oximeter is a high-quality reusable device intended for spot-checking of oxygen saturation and pulse rate in adult patients by using the finger. It is designed for use by both healthcare professionals and lay persons.

We wish you all the best for your health.

Safety notice

Please read these instructions for use carefully before handling the device the first time.

Please store these instructions carefully and ensure other users can access them.


Package contents:

- Pulse oximeter
- 2 x 1.5 V AAA batteries
- Storage bag & lanyard
- Instructions for use


Table of contents	Page
1. Device and display description	88
2. General description	90
3. Measurement principle	91
4. Safety information	92
5. Device features	97
6. Intended use	97
7. Initial operation of the device	98
8. Cleaning and disinfecting	101
9. Maintenance and storage	102
10. Technical data	103
11. Troubleshooting	106
12. Electromagnetic compatibility	108
13. Warranty conditions	108
14. Contact details for customer queries	109
15. Signs and symbols	110

1. Device and display description

1. Symbol of oxygen saturation
2. SpO₂ value
3. Power indication
4. Pulse rate value
5. ON button
6. Perfusion index
7. SpO₂ % Plethysmograph


The pulse plethysmogram displays the user's pulse beat. The height of the plethysmogram shows the user's pulse strength.

The pulse oximeter supports multi-directional display. The screen can display the measured value, pulse wave, and bar graph.

(Figure 1)

When the oxygen saturation and pulse rate signals are incomplete, the measured values are not displayed (Figure 1).

 The bar graph (left part in Figure 1) acts as an indicator of signal incompleteness – the higher the bar graph the stronger the pulse. When the bar graph is not beating smoothly, the accuracy of the measured value may be reduced at this time. When the bar graph is beating smoothly, the measured value read is the best value.

Components

The pulse oximeter consists of probe, electronic circuits, and display and plastic enclosures.

- The probe is the hole in the middle of the equipment to which the finger insert.
- The probe is the applied part of the equipment.

2. General description

Blood Oxygen Saturation (SpO₂)

SpO₂ stands for peripheral capillary oxygen saturation. Oxygen saturation is defined as the ratio of oxyhaemoglobin (HbO₂) to the total concentration of haemoglobin in the blood (i.e. oxyhaemoglobin + reduced haemoglobin). It is an important physiological parameter involved in respiration and circulation. The pulse oximeter is small, portable, non-invasive and easy to use. The user only needs to insert a finger into the chamber to measure the SpO₂ value and pulse rate. The device is intended for use in a home healthcare environment as well as for clinical use.

Reference Values for SpO₂

SpO ₂ measurement range	Classification / Recommendation
95 % – 100 %	Normal Oxygen Level: Indicates sufficient oxygen levels in the blood.
90 % – 94 %	Low Oxygen Level: Keep an eye on your symptoms. A consultation with a physician is recommended.
< 90 %	Critical Oxygen Level: A consultation with a physician is strongly recommended.

Perfusion Index (PI)

Perfusion Index (PI) measures blood flow at the sensor site, indicating the strength of peripheral circulation, i.e., the amount of blood flowing through the body in a given period of time.

It is calculated as the ratio of pulsatile (arterial) blood flow to non-pulsatile (venous and capillary) blood flow.

Reference Values for PI

The PI value range is 0.02 % when the pulse is very low and indicates a very strong pulse at 20 %.

It is important to note that PI is a relative value and can vary significantly based on the individual, measurement site, and environmental conditions.

For this reason, it is recommended that you find out your own normal PI and then use this as a guide when taking measurements. A very low PI value can influence the measurement.

3. Measurement principle

Oxygenated blood absorbs light predominantly at 905 nm (near infrared light), whereas deoxygenated blood absorbs light predominantly at 660 nm (red light). A pulse oximeter passes a beam of red and infrared light through a pulsating capillary bed and then measures the amount of red and infrared light emerging from the tissues via a sensor. To improve accuracy, the pulse oximeter uses a proprietary algorithm to collect data from pulsatile arterial blood and exclude local noise from the tissue. The relative absorption of light by oxyhaemoglobin (HbO₂) and deoxyhaemoglobin is then calculated according to the Beer-Lambert law and a quantitative measurement of the user's oxyhaemoglobin status, i.e. oxygen saturation level (SpO₂), is derived. Due to the sensitivity of the pulse oximeter, the finger should be kept stationary during measurement. It is recommended that you use this device for measurement before or after sports. Do not use for continuous monitoring.

4. Safety information

Caution

(Indicates a potential hazard or unsafe practice that, if not avoided, could result in minor personal injury or device/property damage.)

1. This device does not have an alarm function.
2. Prolonged use or user conditions may require periodic changes to the sensor position. Change the sensor position at least every 2 hours and check skin integrity, circulation condition, and correct alignment.
3. Although the part expected to come into contact with the human body has been biologically evaluated and the biological safety meets the standard requirements, very few people may have allergic reactions. Patients with allergic reactions should stop using the device immediately.
4. The accuracy of pulse oximetry measurements may be adversely affected by high ambient light. Please avoid direct light, such as surgical lights or sunlight, during measurements.
5. Do not disassemble, install, or repair the device without authorisation. Do not look directly at light-emitting devices. It may cause eye damage.
6. The maximum temperature of the blood oxygen probe at the contact surface with the human body does not exceed 41 °C.
7. The data update period is less than 30 seconds. The data update and signal processing affect the display of the SpO₂ value and pulse rate: The blood oxygen (SpO₂) and pulse rate (PR) values displayed on the pulse oximeter are the values obtained by calculating and averaging multiple sets of data acquired within a set period of time. The data is calculated once for each pulse wave. It is a single real-time blood oxygen and pulse rate value. The data averaging method is related to the pulse rate value:
 - when the pulse rate value is below 50 bpm, the last 12 seconds of data are averaged;

- when the pulse rate value is between 50 bpm and 120 bpm, the last 8 seconds of data are averaged;
 - when the pulse rate value is greater than 120 bpm, the last 6 seconds of data are averaged.
- The displayed blood oxygen and pulse rate values are displayed every 1 second.
8. Do not use the device during exercise. The device is not suitable for continuous monitoring.
 9. Do not use for hypoperfusion oximetry.
 10. There are no specific contraindications for pregnant women. Consultation with a physician is recommended.

**Warning**

(Indicates a potential hazard or unsafe practice that, if not avoided, will result in death or serious injury.)

1. Explosive hazard – Do not use the device in the environment with tinder such as anaesthetic.
2. Do not use the device while the patient is being scanned by MRI or CT.
3. The pulse oximeter is NOT intended for neonates and infants.
4. The user's finger thickness should be between 8 to 25.4 mm.
5. It is recommended that the device should be inspected before use. If there is obvious damage, stop using it.

- 
- 
6. Discomfort or pain may occur if the pulse oximeter is used continuously. Especially for patients with a microcirculation barrier, the pulse oximeter should not be used on the same finger more than 10 minutes.
 7. For patients who need a more careful examination of the measurement site, the pulse oximeter should not be placed on oedema or fragile areas.
 8. The pulse oximeter is just a clinical diagnosis auxiliary equipment. The physiological data displayed are for reference only and cannot be directly used for diagnostic interpretation.
 9. It is not recommended to use the device in high frequency environment such as electrosurgical equipment.
 10. Never immerse this device in water or other liquids.
 11. Prevent children from swallowing the device or its accessories. Children must be accompanied by adult guardian when using the device.
 12. Do not service and maintain while the device is in use.
 13. Users are not permitted to maintain the device by themselves.
 14. There are no replaceable components in the device.

**Caution**

(Indicates a potential hazard or unsafe practice that, if not avoided, could result in minor personal injury or device/property damage.)

1. The pulse oximeter is designed to measure the percentage of arterial oxygen saturation of functional haemoglobin. Factors that may degrade the performance of the pulse oximeter or affect the accuracy of the measurement are the following:
 - The pulse oximeter is applied on the same arm as a blood pressure cuff, arterial catheter or infusion line(s) (IVs).
 - Excessive light, such as sunlight or direct home lighting.
 - Moisture in the device.
 - The finger is beyond recommended size range.
 - Poor pulse quality.
 - Venous pulsations.
 - Anaemia or low haemoglobin concentrations.
 - Carboxyhaemoglobin.
 - Methaemoglobin.
 - Dysfunctional haemoglobin.
 - Artificial nails or fingernail polish.
 - Excessive movement during testing may result in incorrect readings.
 - Improper wearing position may result in incorrect readings.
2. Do not look at the light emitted from the device. It is harmful to the eyes.

The patient is an intended operator. All functions can be safely used.

Information about batteries

Power supply (batteries)

- Risk of explosion! Risk of fire! Failure to comply with the following points can result in personal injury or cause overheating, leakage, venting, breakage, explosion, or fire on the battery.
- Batteries must not be charged or short-circuited.
- Do not throw batteries into a fire.
- Never charge, forcibly discharge, heat, disassemble, open, crush, deform, encapsulate, or modify batteries.
- Protect the batteries from direct sunlight, rain, heat, and water.
- Exposure of batteries to an environment with extremely high temperatures or an extremely low air pressure may result in explosion or leakage of flammable liquids and gases.
- Dispose of defective and discharged batteries immediately and properly (see section 'Disposal information').
- Do not use modified or damaged batteries.
- Always select the correct battery type.
- Always insert the batteries correctly, taking into account the polarity (+ / -).
- Never mix old and new batteries or batteries from different manufacturers.
- If a battery has leaked, wear protective gloves and clean the battery compartment with a dry cloth.
- If liquid from a battery cell comes into contact with skin or eyes, clean the affected area with water and seek medical attention if necessary.
- Choking hazard! Keep batteries out of the reach of children. Seek medical attention immediately if swallowed. Swallowing them may cause burns, severe internal injuries, and death.
- Never allow children to replace batteries without adult supervision.

Disposal information

This product is subject to the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment and is marked accordingly. Never dispose of electronic equipment in your household waste. Please obtain information about your local regulations on the proper disposal. This protects the environment and human health.

The batteries comply with the Regulation (EU) 2023/1542.

To protect our environment, empty batteries must not be disposed of in household waste. Please comply with the relevant waste disposal regulations or use public collection points, recycling points or electronic retailers.

5. Device features

- Lightweight, portable and easy to use
- OLED screen shows pulse rate, SpO₂ and perfusion index
- Large font display
- Low battery indicator
- Automatic switch-off function if no signal is detected within 10 seconds
- Sound indicator

6. Intended use

Intended purpose

The Veroval® Pulse oximeter is a reusable device and intended for spot-checking of oxygen saturation and pulse rate. It is to be used on the finger of adults.

Intended user

The Veroval® Pulse oximeter could be used by healthcare professionals and lay persons.

Intended patients

The Veroval® Pulse oximeter is intended for spot-check, non-invasive measurement of arterial oxygen saturation (SpO₂), perfusion index and pulse rate in adult patients. It is suitable for use in both healthy individuals and those with cardiovascular or respiratory conditions, regardless of perfusion status or skin tone.

Indications

The pulse oximeter can effectively estimate and predict the blood oxygen concentration in the periphery of the body by non-invasive operation. With the pulse oximeter changes of respiratory diseases can be monitored by using the finger pulse oxygen as the monitoring index.

Contraindications

The pulse oximeter is not intended for the following cases:


- Do not use the pulse oximeter if there is tissue damage at the measurement site.
 - Do not use to monitor blood oxygen levels in hyperkinetic patients.
 - Do not use to measure blood oxygen levels in hypoperfusion.
- 

7. Initial operation of the device

Battery installation

1. Remove the battery compartment.
2. Insert both batteries into the battery compartment according to the specified polarity.
3. Only use batteries of type: 2 x 1.5 V AAA batteries (LR03). Batteries are included in the delivery scope.
4. Close the battery compartment.

Veroval® Pulse oximeter is now ready for operation.

- If polarity is not correctly aligned, the device will not function and the batteries may leak.
- If the device is not going to be used for a long period of time, remove the batteries from the battery compartment.
- This device should be kept away from children, because small parts can cause children to choke.

Attaching the lanyard

Thread the thin end of the lanyard through the hanging loop on the device, then thread the thick end of the lanyard through the threaded end before tightening it.

Settings

After the pulse oximeter is turned on by pressing the ON button, press and hold the ON button for one second to enter the settings interface.

EN English

Press the ON button once to switch among the settings, press the button for one second to change the current settings; e.g. sounds on or off and other parameters.

Remind Setup	*
Sound Reminder	on
Beep	off
Demo	off
Restore	ok
Brightness	4
Exit	

Limit Setup	*
SpO2 Hi	100
SpO2 Lo	94
PR Hi	130
PR Lo	50
+/-	+
Exit	

1. A beep sounds when the blood oxygen level falls below the set threshold to alert the user.
2. A beep sounds when the pulse rate exceeds the set threshold to alert the user.

The sound function is only used to remind the user to check their blood oxygen level and is not used for any alarm purposes. Any time users have doubts about their blood oxygen saturation, they should check their blood oxygen measurement in time.

Initial use

1. Insert one of your fingers into the finger chamber of the pulse oximeter.
2. Press the ON button to turn on the pulse oximeter.
3. Finger and body should not tremble during measuring.
4. Read the data displayed on the screen.

- When your finger is inserted into the pulse oximeter, your nail surface must be upward.
- The results may be wrong if you do not place your finger correctly in the pulse oximeter.
- Please use medical alcohol to clean the silicon which touches the finger in the pulse oximeter. Also use alcohol to clean the test finger before and after each test. The silicone used in the pulse oximeter is a medical grade silicone. It has no toxins and does not harm the skin.

8. Cleaning and disinfecting

Cleaning

The pulse oximeter should be cleaned on a regular basis. Recommended cleaning agents are ethanol (70 %) or other medical cleanser.

To clean the pulse oximeter, follow these instructions:

- Make sure that the device is turned off.
- Clean the display screen using a soft, clean cloth dampened with a screen cleaner.
- Clean the exterior surface of the device and probe using a soft cloth dampened with the cleaner.
- Wipe off all the cleaning solution with a dry cloth, if necessary.
- Dry the device in a ventilated and cool place.

To avoid damage to the pulse oximeter, follow these instructions:

- Always dilute according to the manufacturer's instructions or use the lowest possible concentration.
- Never immerse this device in water or other liquids.
- Never pour liquid onto the device or accessories.

- Never use abrasive materials, such as steel wool or silver polish, or erosive cleaners, such as acetone or acetone-based cleaners.

Disinfecting

Clean the pulse oximeter before disinfecting it. The recommended disinfectant is ethanol 70 %. The disinfection steps are the same as the cleaning steps. Do not disinfect the device by using high-temperature/high-pressure disinfecting gas.

Never use ethylene oxide or formaldehyde for disinfection.

9. Maintenance and storage

- Replace the batteries in time when the low battery indicator lights up.
- Clean the surface of the pulse oximeter before it is used.
- If the device is not going to be used for a long period of time, remove the batteries from the battery compartment.
- It is best to store the device in a place where ambient temperature is -20 °C to +50 °C and relative humidity is 10 % to 95 %.
- The device should be kept in a dry environment at any time. A wet ambient might affect its lifetime and even might damage the device.
- Avoid exposure to direct sunlight.
- Avoid excessive radioactive rays, infrared rays or ultraviolet rays.

10. Technical data

Model	Veroval® Pulse oximeter	
Type	OXI-70E	
Display	OLED	
SpO ₂	Measurement range	0 – 100 %
	Display range	35 – 100 %
	Accuracy	70 % – 100 % +/-2 %; < 70 % is unspecified
	Resolution	1 %
Pulse rate	Measurement range	25 – 250 bpm
	Display range	25 – 250 bpm
	Accuracy	± 2 bpm
	Resolution	1 bpm
Pulse intensity	Bargraph indicator	
Perfusion index	Measurement range	0 %~50 %
	Display range	0 %~50 %
	Accuracy	0.1 % (0 %~1 %); 1 % (1 %~20 %); > 20 % no definition
	Resolution	1 %
Weight	Approx. 40 g (without batteries)	

Dimension	Approx. 63.7(L) x 36.3(W) x 36.0(H) mm	
Operating conditions	Temperature	+10 °C – +40 °C
	Relative humidity	15 – 95 %, non-condensing
	Atmospheric pressure	700 hPa to 1060 hPa
Storage/transport conditions	Temperature	-20°C – +50°C
	Relative humidity	10 – 95 %, non-condensing
	Atmospheric pressure	700 hPa to 1060 hPa
Peak wave lengths	Red light 660 nm; Infrared 905 nm	
Optical output power	Max. 1.2 mW	
Power supply	2 x 1.5 V alkaline/manganese (AAA/LR03) batteries	
Power consumption	30 mA (normal)	
Battery lifetime	1 year of operation at 4 measurements per day	
Low battery indication	< 1.9 V $\pm$ 0.2 V, automatic switch-off of the device	
Protection against electric shock	Internally powered medical electrical equipment (only when batteries are used) Applied part: Type BF (non-defibrillation proof)	

Protection against harmful ingress of water or solid materials	IP22
Operation mode	Continuous operation
Service life (operating life)	5 years
Reference to standards	The device complies with the EMC requirements of European standard IEC 60601-1-2 (Group 1, Class B, in accordance with CISPR 11, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-8)

i The materials which the user can come into contact with have no toxicity and no action on tissues comply with ISO10993-1, ISO10993-5 and ISO10993-10.

11. Troubleshooting

Error	Possible causes	Remedy
The pulse oximeter cannot be turned on.	The batteries are flat or almost empty.	Check the batteries. If necessary, replace with two identical batteries.
	The batteries are not inserted correctly.	Reinsert the batteries.
	The pulse oximeter is faulty.	Please contact your specialist dealer or customer service.
The display suddenly switches off.	The pulse oximeter is set to switch off automatically in 10 seconds when there are no correct physiological signals.	No action needed.
	The batteries are almost empty.	Replace the batteries.
SpO ₂ or PR are not displayed correctly.	The thoroughfare from photo detector to light emitting diode was sheltered by some objects.	Check and clean the inside of the probe, especially the two sensor windows.
	Finger is trembling or the patient is moving.	Do not move during measurement.

	Finger is not positioned correctly.	Place the finger correctly and try again.
	Finger's size is too big or too small.	Select the correct finger size for the measurement.
	Excessive ambient light	Avoid excessive ambient light irradiation.
	Pulse rate value subject to the cyclical fluctuations.	The patient may have an irregular heart rhythm.
SpO ₂ or PR cannot be displayed.	Finger is not positioned correctly.	Place the finger correctly and try again.
The rate can't be displayed normally.	The measured SpO ₂ is too low to be detected.	Try again. If you are sure that the device is working properly, consult a doctor for a medical advice.
		Make sure nothing is restricting your blood flow.

There are no user-serviceable parts inside the pulse oximeter. The cover should only be removed by qualified service personnel. If you are uncertain about the accuracy of any measurement, check the user's vital signs by alternate means. Also make sure the pulse oximeter is functioning correctly.

Do not spray, pour, or spill any liquid on the pulse oximeter, its accessories, connectors, switches, or openings in the enclosure as this may damage the pulse oximeter.

12. Electromagnetic compatibility

The device is suitable for use in all environments listed in these instructions for use, including domestic environments.

The device complies with the EMC requirements of European standard IEC 60601-1-2 (Group 1, Class B, in accordance with CISPR 11, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-8).

Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. However, if it is necessary to use the device in the manner described above, this device and the other devices must be monitored to ensure that they are working properly.

The use of accessories other than those specified or provided by the manufacturer of this device may result in increased electromagnetic emissions or a decrease in the device's electromagnetic immunity and lead to faulty operation.

Do not use the device near strong electromagnetic fields and keep it away from radio equipment and mobile phones. Portable and mobile high-frequency and communication devices, such as telephones and mobile phones, can impair the functionality of this electronic medical device. Therefore, keep a minimum distance of 30 cm.

Further documentation in accordance with the electromagnetic compatibility standard can be requested from the manufacturer or the customer service (contact details see Chapter 14).

13. Warranty conditions

We offer a 2-year warranty on this premium quality pulse oximeter from the date of purchase and in accordance with the following conditions.

Warranty claims must be made during the warranty period. The purchase date must be documented by an appropriately completed and stamped warranty certificate or proof of purchase.

This device is solely intended for the purpose stated in these instructions for use. The manufacturer is not liable for damage or injuries resulting from improper handling, use or unauthorized interference.

Within the warranty period, HARTMANN will replace or repair free of charge any faulty device components caused by material faults or manufacturing errors. This does not extend the warranty period.

Accessory parts that are subject to wear and tear (batteries, storage bag, packaging etc.) are excluded from warranty. Claims for compensation are limited to the value of the product; compensation for subsequent damage or injury is expressly excluded.

Please send the device, the storage bag, lanyard and the purchase receipt, directly or via your dealer to the Customer Service Department responsible for your country.

14. Contact details for customer queries

ZA – HARTMANN South Africa Epsom Avenue, Northriding,
2169 Johannesburg
www.veroval.info
Tel. +27 860 4278 6266
phzahelpdesk@hartmann.info

If necessary, please contact us at the relevant address listed above for any questions regarding initial startup, use and maintenance of the device, or to report an unexpected operating malfunction or incident.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Signs and symbols

Follow instructions for use

Medical Device

Manufacturer

Distributor

Importer

Date of manufacture

Batch code

Catalogue number

Serial number

Type number

Caution

Safety note indicating possible damage to the device/accessory or a risk of injury or damage to health

Unique Device Identifier

Keep dry

Keep away from sunlight

Temperature limit

Humidity limitation

Atmospheric pressure limitation

Protection against electric shock (type BF)

IP22

Protected against solid objects with a diameter of ≥ 12.5 mm and against vertically falling drops of water when enclosure tilted up to 15°

Direct current

Symbol for the marking of electrical and electronic equipment

European conformity

Dispose of packaging in an environmentally responsible way

Dispose of packaging in an environmentally responsible way

Dispose of packaging in an environmentally responsible way

Authorised representative in the European Community (EC)

Authorised representative in Switzerland

Cardboard recycling code

Alarm inhibit

Do not dispose of batteries containing hazardous substances with household waste

The illustration used in this manual may differ slightly from the appearance of the actual device.

Date of revision of the text: 2025-12-08, Version: A1

Shenzhen AOJ Medical Technology Co., Ltd.
Room 301&4-5F, Block A, Building A, Jingfa Intelligent
Manufacturing Park, Xiawei Yuan,
Gushu Community, Xixiang Street,
Bao'an District, 518126 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA
Email: info@aojmedical.com
Website: <https://www.aojmedical.com>
Tel: 86-755-2778 6026

Share Info GmbH
Am Schulzentrum 12, 41564 Kaarst, Germany
Telephone: 0049 0179 5666 508
Email: EU-Rep@share-info.com

PAUL HARTMANN AG
Paul-Hartmann-Straße 12
89522 HEIDENHEIM, GERMANY
www.hartmann.info

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